

Abking Biotechnologies Inc.
Add: 5F-2, 669, Sec 5, Chunghsiao E. Road, Taipei 110, Taiwan.
Tel: 886-2- 8785-8136
Fax: 886-2- 8785-8103
E-mail: abking@abk.com.tw

Product Information

Kit-2 of FLAXNC-tagged Ab

Cat. Number **ABK-FH001-02**
contains: Peroxidase-Mouse Anti-FLAXNC
monoclonal Ab (Cat. FH0011) and Rabbit
anti-FLAXNC Polyclonal Ab (Cat. Fp0001)

Product Description

Mouse Anti-FLAXNC (FLAG) monoclonal antibody (Cat. Fm0001) is a purified immunoglobulin, IgG1, monoclonal antibody isolated from mouse ascites fluid. Rabbit anti-FLAXNC polyclonal Ab (Cat.Fp0001) is purified from rabbit antiserum. Both Ab react specifically with FLAG tagged recombinant fusion proteins expressed in transfected mammalian cells, from *E. coli* or from in vitro translation. Peroxidase-conjugated Mouse anti-FLAXNC monoclonal antibody is useful in ELISA, Immunoblotting. Rabbit anti-FLAXNC polyclonal antibody is also useful in ELISA, Immunoblotting, Immunofluorescent staining and Immunoprecipitation.

Reagent

Supplied in Freeze-dried powder. Resuspended in 0.4 ml 50% glycerol.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage

Store undiluted antibody at -20°C . Repeated freezing and thawing is not recommended.

Preparation Instructions

Dilute the antibody to 0.2 to $1\ \mu\text{g/ml}$ in 0.05 M Tris buffered saline (TBS), pH 7.4 or PBS. Adjust the antibody concentration to maximize detection sensitivity.

Product Profile

Antigenic binding site: N-Asp-Tyr-Lys-Asp-Asp-Asp-Asp-Lys-C

For monoclonal anti-FLAXNC Ab

Clonality	monoclonal
Clone number	2F11-4H10
Isotype	IgG ₁
Light chain type	kappa
Purity	Protein G purified

Procedures

● Procedure for Western Blot

1. Transfer the FLAG fusion protein of interest from *E. coli* or crude mammalian cell lysate of interest to a PVDF membrane.
2. Block the membrane by using a solution of 5% non-fat dry milk in TBS or PBS at room temperature for 5 minutes to 1 hour.
3. Wash the membrane three times for 5 minutes each in TBS or PBS with 0.05% Tween 20 (Wash buffer-WB) at room temperature.

4. Incubate the membrane with Peroxidase-conjugated Mouse anti-FLAXNC monoclonal antibody (Cat. FH0011) as the primary antibody at 0.2 to $1\ \mu\text{g/ml}$ with 5% non-fat dry milk of Wash Buffer-WB at room temperature for 60 minutes (do not exceed 2 hours). Adjust the antibody concentration to maximize detection sensitivity and to minimize background.

5. Wash the membrane three times for 5 minutes each in Wash buffer-WB at room temperature with agitation.

6. Treat the membrane with chemiluminescent, or other peroxidase substrate.

● Procedure for Immunofluorescent staining

1. Culture mammalian cells on cover slip for proper time. Fix the cover slip with cold acetone/ methanol (1:1) for 90 seconds or 4% paraformaldehyde for 15 minutes.

2. Dilute the antibody to 0.4 to $1\ \mu\text{g/ml}$ in PBS with 0.05% Tween 20 (Wash buffer-IF). After soaking in PBS with 2.5% BSA, cells were stained with rabbit Anti-FLAG polyclonal antibody at 37°C for 30 minutes to overnight.

3. Wash the cell by incubating in Wash buffer-IF at room temperature for three times, each 5 minutes.

4. Incubate the cells with FITC Anti-Rabbit IgG as the secondary antibody in Wash buffer-IF at room temperature for 30 minutes to 1 hour.

5. Wash the cover slip for three times as previous steps. Mount the cells on 50% to 90% glycerol and visualize the cells on fluorescent microscope.

● Procedure for Immuno-precipitation

1. Incubate 0.1 to 1 mg crude cell extract with 2 to $20\ \mu\text{g}$ Rabbit Anti-FLAXNC polyclonal Ab (Cat. Fp0001) in 0.5 ml buffer with proper protease inhibitors at 4°C for 1 hour with agitation.

2. Add $15\ \mu\text{l}$ of 50% slurry of Protein A Sepharose bead and incubate for 30 minutes with agitation.

3. Centrifuge the vial 7000 rpm for 1 minute, and carefully aspirate the supernatant without disturbing the beads.

4. Add 0.5 to 1 ml PBS with 0.1 % Triton (Wash buffer-IP), after inverting for several times or for 3 minutes, centrifuge again as Step 3.

5. Repeat Step 3 & 4 for four to six times.

6. Add sample dye and blot onto PVDF membrane as western blot analysis.

7. After blocking in 5% non-fat dry milk in Wash buffer-IP, use proper primary antibody such as Peroxidase-conjugated Mouse anti- FLAXNC monoclonal antibody (Cat. FH0001) in Wash buffer-IP and incubate for 45 minutes.

8. Treat the membrane with chemiluminescent, or other peroxidase substrate.